



*In freezing temperatures, searing heat and even underwater...you can rely on NJ Phillips.*



## Draw-Off Caps and Spikes

Code: see below



**Draw-Off**  
(Specific to Merial)

- **WX230** – Designed to suit 20mm bottle neck (pack of 2).
- **WX231** – Designed to suit 30mm bottle neck (pack of 2).
- **WX1345** – Designed to suit 32mm bottle neck (pack of 2).
- **WX1501** – Designed for use with Virbamec and Virbamec Plus.



**Draw-Off**  
Code: WX1501  
(Specific to Virbac)

## Vented Draw-Off Cap

Code: WX1319A



- Developed for Virbac USA to be used with Virbamec Pour-On.
- Dual-barb spigot 1/4" to 3/8".
- Vented to allow air into bottle for ease of product application and applicator efficiency.

SPECIFIC TO VIRBAC

## Draw-Off Spike

Code: WX1637

- 160mm draw-off spike to suit large bottles (pack of 1).

